

Region/View UPPER EXTREMITY		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High	
L R		0		Medium		0		0	
Apr Ref		Patient		kV	mA	Sec	mAs	Index	Cm
1.1.4		Small		80	350	.090	31.5	4	4
1.1.5		Medium		80	350	.090	31.5	5	5
1.1.6		Large		80	350	.090	31.5	6	6

Region/View UPPER EXTREMITY			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High		
L R		0		Medium		0		0		
Apr Ref		Patient		kV	mA	Sec	mAs	Index		Cm
1.2.4		Small		80	350	.090	31.5	4		4
1.2.5		Medium		80	350	.090	31.5	5		5
1.2.6		Large		80	350	.090	31.5	6		6

Region/View UPPER EXTREMITY			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High		
L R		0		Medium		0		0		
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
1.3.4		Small		80	350	.090	31.5		4	4
1.3.5		Medium		80	350	.090	31.5		5	5
1.3.6		Large		80	350	.090	31.5		6	6

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
UPPER EXTREMITY								
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
1.4.4	Small	80	350	.090	31.5	4	4	
1.4.5	Medium	80	350	.090	31.5	5	5	
1.4.6	Large	80	350	.090	31.5	6	6	

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
UPPER EXTREMITY							
FINGER		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
1.5.4	Small	52	100	.030	3	4	4
1.5.5	Medium	53	100	.030	3	5	5
1.5.6	Large	53	100	.035	3.5	6	6

Region/View UPPER EXTREMITY				Focal Spot		Receptor		AEC Enabled		Film Screen		
HAND				S		None		N		2		
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High			
L		-50		Medium		0			0			
Apr Ref		Patient		kV	mA	Sec		mAs		Index		Cm
1.6.4		Small		53	100	.050		5		4		4
1.6.5		Medium		54	100	.050		5		5		5
1.6.6		Large		54	100	.060		6		6		6

Region/View		Focal	Receptor	AEC	Film		
UPPER	EXTREMITY	Spot		Enabled	Screen		
WRIST		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low	Cm Limit High			
L	-50	Medium	0	0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
1.7.4	Small	51	100	.060	6	4	4
1.7.5	Medium	51	100	.080	8	5	5
1.7.6	Large	52	100	.080	8	6	6

Region/View		Focal		Receptor		AEC		Film	
UPPER EXTREMITY		Spot				Enabled		Screen	
FOREARM		S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High	
L		-50		Medium		0		0	
Apr Ref		Patient		kV mA		Sec mAs		Index Cm	
1.8.5		Small		57 100		.060 6		5 5	
1.8.6		Medium		58 100		.080 8		6 6	
1.8.8		Large		60 100		.090 9		8 8	

Region/View UPPER EXTREMITY		Focal Spot		Receptor		AEC Enabled		Film Screen	
ELBOW		S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
1.9.5	Small	57	100	.090	9	5	7		
1.9.6	Medium	59	100	.090	9	6	8		
1.9.7	Large	59	100	.100	10	7	9		

Region/View UPPER EXTREMITY		Focal Spot		Receptor		AEC Enabled		Film Screen	
HUMERUS		S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High	
L		-50		Medium		0		0	
Apr Ref		Patient		kV mA		Sec mAs		Index Cm	
1.10.6		Small		56 100		.100 10		6 8	
1.10.7		Medium		58 100		.100 10		7 9	
1.10.9		Large		60 200		.060 12		9 11	

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
UPPER EXTREMITY							
SHOULDER		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
1.11.6	Small	67	100	.180	18	6	9
1.11.7	Medium	67	100	.180	18	7	10
1.11.9	Large	69	200	.120	24	9	12

Region/View		Focal Spot		Receptor		AEC Enabled	Film Screen
UPPER EXTREMITY							
SPARE		S		None		N	1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
1.12.4	Small	80	350	.090	31.5	4	4
1.12.5	Medium	80	350	.090	31.5	5	5
1.12.6	Large	80	350	.090	31.5	6	6

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
LOWER EXTREMITY								
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
2.1.4	Small	80	350	.090	31.5	4	4	
2.1.5	Medium	80	350	.090	31.5	5	5	
2.1.6	Large	80	350	.090	31.5	6	6	

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
LOWER EXTREMITY								
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
2.2.4	Small	80	350	.090	31.5	4	4	
2.2.5	Medium	80	350	.090	31.5	5	5	
2.2.6	Large	80	350	.090	31.5	6	6	

Region/View LOWER EXTREMITY				Focal Spot		Receptor		AEC Enabled		Film Screen		
SPARE				S		None		N		1		
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High			
L R		0		Medium		0			0			
Apr Ref		Patient		kV	mA	Sec		mAs		Index		Cm
2.3.4		Small		80	350	.090		31.5		4		4
2.3.5		Medium		80	350	.090		31.5		5		5
2.3.6		Large		80	350	.090		31.5		6		6

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen	
LOWER EXTREMITY									
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
2.4.4	Small	80	350	.090	31.5	4	4		
2.4.5	Medium	80	350	.090	31.5	5	5		
2.4.6	Large	80	350	.090	31.5	6	6		

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
LOWER EXTREMITY								
FINGER		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High	
L	-50	Medium		0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
2.5.4	Small	52	100	.030	3	4	4	
2.5.5	Medium	53	100	.030	3	5	5	
2.5.6	Large	53	100	.035	3.5	6	6	

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
LOWER EXTREMITY								
HAND		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L	-50	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
2.6.4	Small	53	100	.050	5	4	4	
2.6.5	Medium	54	100	.050	5	5	5	
2.6.6	Large	54	100	.060	6	6	6	

Region/View		Focal Spot		Receptor		AEC Enabled	Film Screen
LOWER EXTREMITY							
WRIST		S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
2.7.4	Small	51	100	.060	6	4	4
2.7.5	Medium	51	100	.080	8	5	5
2.7.6	Large	52	100	.080	8	6	6

Region/View		Focal	Receptor	AEC	Film		
LOWER	EXTREMITY	Spot		Enabled	Screen		
FOREARM		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
2.8.5	Small	57	100	.060	6	5	5
2.8.6	Medium	58	100	.080	8	6	6
2.8.8	Large	60	100	.090	9	8	8

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
LOWER EXTREMITY							
ELBOW		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low	Cm Limit High			
L	-50	Medium	0	0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
2.9.5	Small	57	100	.090	9	5	7
2.9.6	Medium	59	100	.090	9	6	8
2.9.7	Large	59	100	.100	10	7	9

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
LOWER EXTREMITY							
HUMERUS		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
2.10.6	Small	56	100	.100	10	6	8
2.10.7	Medium	58	100	.100	10	7	9
2.10.9	Large	60	200	.060	12	9	11

Region/View		Focal	Receptor	AEC	Film		
LOWER	EXTREMITY	Spot		Enabled	Screen		
SHOULDER		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low	Cm Limit High			
L	-50	Medium	0	0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
2.11.6	Small	67	100	.180	18	6	9
2.11.7	Medium	67	100	.180	18	7	10
2.11.9	Large	69	200	.120	24	9	12

Region/View LOWER EXTREMITY				Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE				S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High		
L		-50		Medium		0			0		
Apr Ref		Patient		kV	mA	Sec	mAs		Index		Cm
2.12.4		Small		80	350	.090	31.5		4		4
2.12.5		Medium		80	350	.090	31.5		5		5
2.12.6		Large		80	350	.090	31.5		6		6

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
3.1.4	Small	80	350	.090	31.5	4	4		
3.1.5	Medium	80	350	.090	31.5	5	5		
3.1.6	Large	80	350	.090	31.5	6	6		

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L R	0	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs		Index	Cm	
3.2.4	Small	80	350	.090	31.5		4	4	
3.2.5	Medium	80	350	.090	31.5		5	5	
3.2.6	Large	80	350	.090	31.5		6	6	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High	
L R		0		Medium		0		0	
Apr Ref		Patient		kV	mA	Sec	mAs	Index	Cm
3.3.4		Small		80	350	.090	31.5	4	4
3.3.5		Medium		80	350	.090	31.5	5	5
3.3.6		Large		80	350	.090	31.5	6	6

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
3.4.4	Small	80	350	.090	31.5	4	4		
3.4.5	Medium	80	350	.090	31.5	5	5		
3.4.6	Large	80	350	.090	31.5	6	6		

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
FINGER		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs		Index	Cm	
3.5.4	Small	52	100	.030	3		4	4	
3.5.5	Medium	53	100	.030	3		5	5	
3.5.6	Large	53	100	.035	3.5		6	6	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
HAND		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs		Index	Cm	
3.6.4	Small	53	100	.050	5	4		4	
3.6.5	Medium	54	100	.050	5	5		5	
3.6.6	Large	54	100	.060	6	6		6	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
WRIST		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L	-50	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
3.7.4	Small	51	100	.060	6	4	4		
3.7.5	Medium	51	100	.080	8	5	5		
3.7.6	Large	52	100	.080	8	6	6		

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
FOREARM		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L	-50	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
3.8.5	Small	57	100	.060	6	5		5	
3.8.6	Medium	58	100	.080	8	6		6	
3.8.8	Large	60	100	.090	9	8		8	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
ELBOW		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L	-50	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
3.9.5	Small	57	100	.090	9	5		7	
3.9.6	Medium	59	100	.090	9	6		8	
3.9.7	Large	59	100	.100	10	7		9	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled	Film Screen
HUMERUS		S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
3.10.6	Small	56	100	.100	10	6	8
3.10.7	Medium	58	100	.100	10	7	9
3.10.9	Large	60	200	.060	12	9	11

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled		Film Screen	
SHOULDER		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs		Index	Cm	
3.11.6	Small	67	100	.180	18		6	9	
3.11.7	Medium	67	100	.180	18		7	10	
3.11.9	Large	69	200	.120	24		9	12	

Region/View PELVIS & HIP		Focal Spot		Receptor		AEC Enabled	Film Screen
SPARE		S		None		N	1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
3.12.4	Small	80	350	.090	31.5	4	4
3.12.5	Medium	80	350	.090	31.5	5	5
3.12.6	Large	80	350	.090	31.5	6	6

Region/View ABDOMEN K U B		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
4.1.4	Small	80	350	.090	31.5	4		4	
4.1.5	Medium	80	350	.090	31.5	5		5	
4.1.6	Large	80	350	.090	31.5	6		6	

Region/View ABDOMEN K U B		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
4.2.4	Small	80	350	.090	31.5	4		4	
4.2.5	Medium	80	350	.090	31.5	5		5	
4.2.6	Large	80	350	.090	31.5	6		6	

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L R		0		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.3.4		Small		80	350	.090	31.5		4	4
4.3.5		Medium		80	350	.090	31.5		5	5
4.3.6		Large		80	350	.090	31.5		6	6

Region/View ABDOMEN K U B		Focal Spot	Receptor		AEC Enabled		Film Screen
SPARE		S	None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L R	0	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
4.4.4	Small	80	350	.090	31.5	4	4
4.4.5	Medium	80	350	.090	31.5	5	5
4.4.6	Large	80	350	.090	31.5	6	6

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
ABDOMEN K U B										
FINGER			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index		Cm
4.5.4		Small	52	100	.030	3		4		4
4.5.5		Medium	53	100	.030	3		5		5
4.5.6		Large	53	100	.035	3.5		6		6

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
HAND			S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.6.4		Small		53	100	.050	5		4	4
4.6.5		Medium		54	100	.050	5		5	5
4.6.6		Large		54	100	.060	6		6	6

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
WRIST			S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.7.4		Small		51	100	.060	6	4	4	
4.7.5		Medium		51	100	.080	8	5	5	
4.7.6		Large		52	100	.080	8	6	6	

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
FOREARM			S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.8.5		Small		57	100	.060	6	5	5	
4.8.6		Medium		58	100	.080	8	6	6	
4.8.8		Large		60	100	.090	9	8	8	

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
ELBOW			S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.9.5		Small		57	100	.090	9		5	7
4.9.6		Medium		59	100	.090	9		6	8
4.9.7		Large		59	100	.100	10		7	9

Region/View ABDOMEN K U B				Focal Spot		Receptor		AEC Enabled		Film Screen	
HUMERUS				S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High		
L		-50		Medium		0			0		
Apr Ref		Patient		kV	mA	Sec	mAs		Index		Cm
4.10.6		Small		56	100	.100	10		6		8
4.10.7		Medium		58	100	.100	10		7		9
4.10.9		Large		60	200	.060	12		9		11

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
ABDOMEN K U B										
SHOULDER			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm			
4.11.6	Small	67	100	.180	18	6	9			
4.11.7	Medium	67	100	.180	18	7	10			
4.11.9	Large	69	200	.120	24	9	12			

Region/View ABDOMEN K U B			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
4.12.4		Small		80	350	.090	31.5		4	4
4.12.5		Medium		80	350	.090	31.5		5	5
4.12.6		Large		80	350	.090	31.5		6	6

Region/View CERVICAL SPINE		Focal Spot	Receptor	AEC Enabled	Film Screen		
SPARE		S	None	N	1		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L R	0	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
5.1.4	Small	80	350	.090	31.5	4	4
5.1.5	Medium	80	350	.090	31.5	5	5
5.1.6	Large	80	350	.090	31.5	6	6

Region/View CERVICAL SPINE		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
5.2.4	Small	80	350	.090	31.5	4		4	
5.2.5	Medium	80	350	.090	31.5	5		5	
5.2.6	Large	80	350	.090	31.5	6		6	

Region/View CERVICAL SPINE		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
5.3.4	Small	80	350	.090	31.5	4		4	
5.3.5	Medium	80	350	.090	31.5	5		5	
5.3.6	Large	80	350	.090	31.5	6		6	

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
CERVICAL SPINE								
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
5.4.4	Small	80	350	.090	31.5	4	4	
5.4.5	Medium	80	350	.090	31.5	5	5	
5.4.6	Large	80	350	.090	31.5	6	6	

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
CERVICAL SPINE										
FINGER			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index	Cm	
5.5.4		Small	52	100	.030	3		4	4	
5.5.5		Medium	53	100	.030	3		5	5	
5.5.6		Large	53	100	.035	3.5		6	6	

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
CERVICAL SPINE								
HAND		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L	-50	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
5.6.4	Small	53	100	.050	5	4	4	
5.6.5	Medium	54	100	.050	5	5	5	
5.6.6	Large	54	100	.060	6	6	6	

Region/View		Focal Spot		Receptor		AEC Enabled	Film Screen
CERVICAL SPINE							
WRIST		S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
5.7.4	Small	51	100	.060	6	4	4
5.7.5	Medium	51	100	.080	8	5	5
5.7.6	Large	52	100	.080	8	6	6

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen
CERVICAL SPINE								
FOREARM		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High	
L	-50	Medium		0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
5.8.5	Small	57	100	.060	6	5	5	
5.8.6	Medium	58	100	.080	8	6	6	
5.8.8	Large	60	100	.090	9	8	8	

Region/View		Focal Spot		Receptor		AEC Enabled		Film Screen	
CERVICAL SPINE									
ELBOW		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
5.9.5	Small	57	100	.090	9	5		7	
5.9.6	Medium	59	100	.090	9	6		8	
5.9.7	Large	59	100	.100	10	7		9	

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
CERVICAL SPINE							
HUMERUS		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
5.10.6	Small	56	100	.100	10	6	8
5.10.7	Medium	58	100	.100	10	7	9
5.10.9	Large	60	200	.060	12	9	11

Region/View		Focal Spot		Receptor		AEC Enabled	Film Screen
CERVICAL SPINE							
SHOULDER		S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
5.11.6	Small	67	100	.180	18	6	9
5.11.7	Medium	67	100	.180	18	7	10
5.11.9	Large	69	200	.120	24	9	12

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
CERVICAL SPINE										
SPARE			S		None		N		1	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs	Index		Cm	
5.12.4		Small	80	350	.090	31.5	4		4	
5.12.5		Medium	80	350	.090	31.5	5		5	
5.12.6		Large	80	350	.090	31.5	6		6	

Region/View LUMBAR SPINE		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
6.1.4	Small	80	350	.090	31.5	4	4		
6.1.5	Medium	80	350	.090	31.5	5	5		
6.1.6	Large	80	350	.090	31.5	6	6		

Region/View LUMBAR SPINE		Focal Spot		Receptor		AEC Enabled		Film Screen
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
6.2.4	Small	80	350	.090	31.5	4	4	
6.2.5	Medium	80	350	.090	31.5	5	5	
6.2.6	Large	80	350	.090	31.5	6	6	

Region/View LUMBAR SPINE			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields		Density		Patient Size		Cm Limit Low		Cm Limit High		
L R		0		Medium		0		0		
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
6.3.4		Small		80	350	.090	31.5		4	4
6.3.5		Medium		80	350	.090	31.5		5	5
6.3.6		Large		80	350	.090	31.5		6	6

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
LUMBAR SPINE										
SPARE			S		None		N		1	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L R		0	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index		Cm
6.4.4		Small	80	350	.090	31.5		4		4
6.4.5		Medium	80	350	.090	31.5		5		5
6.4.6		Large	80	350	.090	31.5		6		6

Region/View		Focal		Receptor		AEC		Film
LUMBAR	SPINE	Spot				Enabled		Screen
FINGER		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L	-50	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
6.5.4	Small	52	100	.030	3	4	4	
6.5.5	Medium	53	100	.030	3	5	5	
6.5.6	Large	53	100	.035	3.5	6	6	

Region/View LUMBAR SPINE		Focal Spot		Receptor		AEC Enabled		Film Screen
HAND		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High	
L	-50	Medium		0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
6.6.4	Small	53	100	.050	5	4	4	
6.6.5	Medium	54	100	.050	5	5	5	
6.6.6	Large	54	100	.060	6	6	6	

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
LUMBAR	SPINE						
WRIST		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
6.7.4	Small	51	100	.060	6	4	4
6.7.5	Medium	51	100	.080	8	5	5
6.7.6	Large	52	100	.080	8	6	6

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
LUMBAR SPINE										
FOREARM			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index		Cm
6.8.5		Small	57	100	.060	6		5		5
6.8.6		Medium	58	100	.080	8		6		6
6.8.8		Large	60	100	.090	9		8		8

Region/View		Focal		Receptor		AEC		Film
LUMBAR	SPINE	Spot				Enabled		Screen
ELBOW		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L	-50	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
6.9.5	Small	57	100	.090	9	5	7	
6.9.6	Medium	59	100	.090	9	6	8	
6.9.7	Large	59	100	.100	10	7	9	

Region/View		Focal Spot	Receptor	AEC Enabled	Film Screen		
LUMBAR	SPINE						
HUMERUS		S	None	N	2		
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
6.10.6	Small	56	100	.100	10	6	8
6.10.7	Medium	58	100	.100	10	7	9
6.10.9	Large	60	200	.060	12	9	11

Region/View			Focal Spot		Receptor		AEC Enabled		Film Screen	
LUMBAR SPINE										
SHOULDER			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index		Cm
6.11.6		Small	67	100	.180	18		6		9
6.11.7		Medium	67	100	.180	18		7		10
6.11.9		Large	69	200	.120	24		9		12

Region/View LUMBAR SPINE		Focal Spot	Receptor		AEC Enabled	Film Screen	
SPARE		S	None		N	1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
6.12.4	Small	80	350	.090	31.5	4	4
6.12.5	Medium	80	350	.090	31.5	5	5
6.12.6	Large	80	350	.090	31.5	6	6

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled	Film Screen
SPARE		S		None		N	1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L R	0	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.1.4	Small	80	350	.090	31.5	4	4
7.1.5	Medium	80	350	.090	31.5	5	5
7.1.6	Large	80	350	.090	31.5	6	6

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled	Film Screen
SPARE		S		None		N	1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L R	0	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.2.4	Small	80	350	.090	31.5	4	4
7.2.5	Medium	80	350	.090	31.5	5	5
7.2.6	Large	80	350	.090	31.5	6	6

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled	Film Screen
SPARE		S		None		N	1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L R	0	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.3.4	Small	80	350	.090	31.5	4	4
7.3.5	Medium	80	350	.090	31.5	5	5
7.3.6	Large	80	350	.090	31.5	6	6

Region/View CHEST & THORAX		Focal Spot	Receptor		AEC Enabled	Film Screen	
SPARE		S	None		N	1	
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L R	0	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.4.4	Small	80	350	.090	31.5	4	4
7.4.5	Medium	80	350	.090	31.5	5	5
7.4.6	Large	80	350	.090	31.5	6	6

Region/View CHEST & THORAX		Focal Spot	Receptor		AEC Enabled	Film Screen	
FINGER		S	None		N	2	
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.5.4	Small	52	100	.030	3	4	4
7.5.5	Medium	53	100	.030	3	5	5
7.5.6	Large	53	100	.035	3.5	6	6

Region/View CHEST & THORAX		Focal Spot	Receptor		AEC Enabled	Film Screen	
HAND		S	None		N	2	
AEC Fields	Density	Patient Size	Cm Limit Low		Cm Limit High		
L	-50	Medium	0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.6.4	Small	53	100	.050	5	4	4
7.6.5	Medium	54	100	.050	5	5	5
7.6.6	Large	54	100	.060	6	6	6

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled		Film Screen	
WRIST		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
7.7.4	Small	51	100	.060	6	4		4	
7.7.5	Medium	51	100	.080	8	5		5	
7.7.6	Large	52	100	.080	8	6		6	

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled		Film Screen	
FOREARM		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
7.8.5	Small	57	100	.060	6	5		5	
7.8.6	Medium	58	100	.080	8	6		6	
7.8.8	Large	60	100	.090	9	8		8	

Region/View CHEST & THORAX			Focal Spot		Receptor		AEC Enabled		Film Screen	
ELBOW			S		None		N		2	
AEC Fields	Density		Patient Size		Cm Limit Low			Cm Limit High		
L	-50		Medium		0			0		
Apr Ref	Patient		kV	mA	Sec	mAs		Index		Cm
7.9.5	Small		57	100	.090	9		5		7
7.9.6	Medium		59	100	.090	9		6		8
7.9.7	Large		59	100	.100	10		7		9

Region/View CHEST & THORAX		Focal Spot	Receptor		AEC Enabled	Film Screen	
HUMERUS		S	None		N	2	
AEC Fields	Density	Patient Size	Cm Limit Low			Cm Limit High	
L	-50	Medium	0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.10.6	Small	56	100	.100	10	6	8
7.10.7	Medium	58	100	.100	10	7	9
7.10.9	Large	60	200	.060	12	9	11

Region/View CHEST & THORAX		Focal Spot		Receptor		AEC Enabled	Film Screen
SHOULDER		S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.11.6	Small	67	100	.180	18	6	9
7.11.7	Medium	67	100	.180	18	7	10
7.11.9	Large	69	200	.120	24	9	12

Region/View CHEST & THORAX		Focal Spot	Receptor		AEC Enabled	Film Screen	
SPARE		S	None		N	1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High	
L	-50	Medium		0		0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm
7.12.4	Small	80	350	.090	31.5	4	4
7.12.5	Medium	80	350	.090	31.5	5	5
7.12.6	Large	80	350	.090	31.5	6	6

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE		S		None		N		1	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L R	0	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm		
8.1.4	Small	80	350	.090	31.5	4	4		
8.1.5	Medium	80	350	.090	31.5	5	5		
8.1.6	Large	80	350	.090	31.5	6	6		

Region/View SKULL & FACIAL			Focal Spot		Receptor		AEC Enabled		Film Screen	
SPARE			S		None		N		1	
AEC Fields	Density		Patient Size		Cm Limit Low			Cm Limit High		
L R	0		Medium		0			0		
Apr Ref	Patient		kV	mA	Sec	mAs		Index	Cm	
8.2.4	Small		80	350	.090	31.5		4	4	
8.2.5	Medium		80	350	.090	31.5		5	5	
8.2.6	Large		80	350	.090	31.5		6	6	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled	Film Screen	
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
8.3.4	Small	80	350	.090	31.5	4	4	
8.3.5	Medium	80	350	.090	31.5	5	5	
8.3.6	Large	80	350	.090	31.5	6	6	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L R	0	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
8.4.4	Small	80	350	.090	31.5	4	4	
8.4.5	Medium	80	350	.090	31.5	5	5	
8.4.6	Large	80	350	.090	31.5	6	6	

Region/View SKULL & FACIAL			Focal Spot		Receptor		AEC Enabled		Film Screen	
FINGER			S		None		N		2	
AEC Fields		Density	Patient Size		Cm Limit Low			Cm Limit High		
L		-50	Medium		0			0		
Apr Ref		Patient	kV	mA	Sec	mAs		Index		Cm
8.5.4		Small	52	100	.030	3		4		4
8.5.5		Medium	53	100	.030	3		5		5
8.5.6		Large	53	100	.035	3.5		6		6

Region/View SKULL & FACIAL			Focal Spot		Receptor		AEC Enabled		Film Screen	
HAND			S		None		N		2	
AEC Fields	Density		Patient Size		Cm Limit Low			Cm Limit High		
L	-50		Medium		0			0		
Apr Ref	Patient		kV	mA	Sec	mAs		Index	Cm	
8.6.4	Small		53	100	.050	5		4	4	
8.6.5	Medium		54	100	.050	5		5	5	
8.6.6	Large		54	100	.060	6		6	6	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen	
WRIST		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High		
L	-50	Medium		0			0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
8.7.4	Small	51	100	.060	6	4		4	
8.7.5	Medium	51	100	.080	8	5		5	
8.7.6	Large	52	100	.080	8	6		6	

Region/View SKULL & FACIAL			Focal Spot		Receptor		AEC Enabled		Film Screen	
FOREARM			S		None		N		2	
AEC Fields		Density		Patient Size		Cm Limit Low			Cm Limit High	
L		-50		Medium		0			0	
Apr Ref		Patient		kV	mA	Sec	mAs		Index	Cm
8.8.5		Small		57	100	.060	6		5	5
8.8.6		Medium		58	100	.080	8		6	6
8.8.8		Large		60	100	.090	9		8	8

Region/View SKULL & FACIAL			Focal Spot		Receptor		AEC Enabled	Film Screen
ELBOW			S		None		N	2
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High	
L	-50	Medium		0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
8.9.5	Small	57	100	.090	9	5	7	
8.9.6	Medium	59	100	.090	9	6	8	
8.9.7	Large	59	100	.100	10	7	9	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen
HUMERUS		S		None		N		2
AEC Fields	Density	Patient Size		Cm Limit Low			Cm Limit High	
L	-50	Medium		0			0	
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
8.10.6	Small	56	100	.100	10	6	8	
8.10.7	Medium	58	100	.100	10	7	9	
8.10.9	Large	60	200	.060	12	9	11	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen	
SHOULDER		S		None		N		2	
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High			
L	-50	Medium		0		0			
Apr Ref	Patient	kV	mA	Sec	mAs	Index		Cm	
8.11.6	Small	67	100	.180	18	6		9	
8.11.7	Medium	67	100	.180	18	7		10	
8.11.9	Large	69	200	.120	24	9		12	

Region/View SKULL & FACIAL		Focal Spot		Receptor		AEC Enabled		Film Screen
SPARE		S		None		N		1
AEC Fields	Density	Patient Size		Cm Limit Low		Cm Limit High		
L	-50	Medium		0		0		
Apr Ref	Patient	kV	mA	Sec	mAs	Index	Cm	
8.12.4	Small	80	350	.090	31.5	4	4	
8.12.5	Medium	80	350	.090	31.5	5	5	
8.12.6	Large	80	350	.090	31.5	6	6	