

# MELTI Piston



- Large Loop For Easy Placement On Incus
- 2x Wide Standard Platinum Ribbon Prevents Rolling When Crimping
- 0.5mm Depth Marker To Aid Surgeon In Visualizing Placement Into the Vestible
- Smooth Transitioned Depth Marker Prevents Interference With Fenestration In Footplate

Cat. #	Materials	Shaft Dia.	Length
418-375	Platinum & Titanium	0.5mm	3.75mm
418-400	Platinum & Titanium	0.5mm	4.00mm
418-425	Platinum & Titanium	0.5mm	4.25mm
418-450	Platinum & Titanium	0.5mm	4.50mm
418-475	Platinum & Titanium	0.5mm	4.75mm

Designed in conjunction with Dennis Bojrab, MD and Michael LaRouere, MD, Detroit, MI