

ALTO NANO

Need a little more room...downsize to an accomodating design...the ALTO Nano

Partial Ref.#653



Total Ref.#603



Nano platform
2.5mm x 1.5mm



Traditional platform
3.5mm x 2.5mm



- **Downsized** platform allows for exact placement despite surrounding bone.
- Second look **benefits** from the mastoid...reduced need to lift the TM flap.
- **Fits** nicely between malleus & EAC if scutum is intact.
- Smaller size allows covering only the scutal area **reducing chances** of extrusion around the cartilage graft.
- Smaller size **supports** cartilage graft with enhanced contact.
- Packaged in 'Adjuster System' for easy sizing.



New! mm scale for **EASY** measurement — 3 and 6mm engraved numerals for **QUICK** reference



*Designed in conjunction with Peter Weber, MD, Cleveland, OH

GRACE MEDICAL

an enteroptyx company

8500 Wolf Lake Dr., Suite 111 • Memphis, TN 38133 USA • 1.866.472.2363 • 901.386.0990 • Fax 901.386.0950

www.gracemedical.com • info@gracemedical.com